



Tel: + 358 (0)207 917 600
Fax: + 358 (0)207 917 649
E-mail: repro.jt@jt.fi
www.jt.fi

JAAKKO-TAARA OY
RUISSALONTIE 4
20200 TURKU
FINLAND

Colours are referential.

NOTE! We will charge changes
that are against manuscript.

Please add changes
as clearly as possible.

JT is not responsible of
mistakes in approved proof.

JT Job No.: 258006

Item No.: 81911002

Drawn up/OS: Marianne Murto

Date: 26.9.2013

Client: Bayer Oy

Dimensions (mm):

50 x 120

Reeling: 2

Font / Font size / Spacing:

BayerSans / 6,5 pt / 7,15 pt

Colors: Pms 321, warm red, black

Uncoated area: -

Pharmacode: 2377



Approved for printing / Date



Mirena®
Levonorgestrel

20 micrograms/24 hours
intrauterine delivery system

1 sterile intrauterine delivery
system contains:
Levonorgestrel. 52 mg
Constit. q.s.

Store below 30°C
**Keep out of the reach and sight of
children**

To be inserted by a health care professional
Sterile unless the blister is damaged or open

- ▶ **For loading, push the slider
upwards**
- ▶ **This system cannot be re-loaded**
- ▶ **Please read the complete insertion
instructions before insertion**

Manufacturer:
Bayer Oy
Pansiontie 47,
20210 Turku, Finland

Harus dengan resep dokter

Reg. No.: DK10775601476A1

MAL19986422AZ
Controlled Medicine/Ubat Terkawal
SIN 07915P
BRU 11070730P



81911002

Lot No:

Mfg Dt:

To be inserted before:
(อย่าลืมอายุ)



Insertion Instructions

To be inserted by a health care professional using aseptic technique.

Mirena is supplied within an inserter in a sterile package which should not be opened until needed for insertion. Do not resterilize. As supplied, Mirena is for single use only. Do not use if the inner package is damaged or open. Do not insert after the expiry month and year shown on the label.


For timing of insertion, please consult the Mirena prescribing information.

Preparation for insertion

- ▶ Examine the patient to establish the size and position of the uterus, in order to detect any signs of acute genital infections or other contraindications for the insertion of Mirena and to exclude pregnancy.
- ▶ Insert a speculum, visualize the cervix, and then thoroughly cleanse the cervix and vagina with a suitable antiseptic solution.
- ▶ Employ an assistant as necessary.
- ▶ Grasp the anterior lip of the cervix with a tenaculum or other forceps to stabilize the uterus. If the uterus is retroverted, it may be more appropriate to grasp the posterior lip of the cervix. Gentle traction on the forceps can be applied to straighten the cervical canal. The forceps should remain in position and gentle counter traction on the cervix should be maintained throughout the insertion procedure.
- ▶ Advance a uterine sound through the cervical canal to the fundus to measure the depth and confirm the direction of the uterine cavity and to exclude any evidence of intrauterine abnormalities (e.g., septum, submucous fibroids) or a previously inserted intrauterine contraceptive which has not been removed. If difficulty is encountered, consider dilatation of the canal. If cervical dilatation is required, consider using analgesics and/or a paracervical block.

Insertion


- 1.** First, open the sterile package completely (Figure 1). Then use sterile technique and sterile gloves.



- 2.** Push the slider **forward** in the direction of the arrow to the furthest position to load Mirena into the insertion tube (Figure 2).

IMPORTANT! Do not pull the slider downwards as this may prematurely release Mirena. Once released, Mirena cannot be re-loaded.

- 3.** Holding the slider in the furthest position, set the **upper** edge of the flange to correspond to the sound measurement of the uterine depth (Figure 3).






- 4.** While holding the slider in the **furthest** position, advance the inserter through the cervix until the flange is approx. 1.5-2.0 cm from the uterine cervix (Figure 4).

IMPORTANT! Do not force the inserter. Dilate the cervical canal, if necessary.

- 5.** While holding the inserter steady, **pull the slider to the mark** to open the horizontal arms of Mirena (Figure 5). Wait 5 -10 seconds for the horizontal arms to open completely.

- 6.** Advance the inserter gently towards the fundus of the uterus **until the flange touches the cervix**. Mirena is now in the fundal position (Figure 6).



- 7.** Holding the inserter in place, release Mirena by pulling **the slider all the way down** (Figure 7). While holding the slider all the way down, gently remove the inserter by pulling it out. **Cut the threads** to leave about 2-3 cm visible outside of the cervix.


IMPORTANT! Should you suspect that the system is not in the correct position, check placement (e.g. with ultrasound). Remove the system if it is not positioned properly within the uterine cavity. A removed system must not be re-inserted.

Removal / replacement

For removal/replacement, please consult the Mirena prescribing information.

Mirena is removed by pulling on the threads with a forceps (Figure 8).

You may insert a new Mirena immediately following removal.



MANUFACTURED BY:

Bayer Oy
Pansiontie 47
20210 Turku
Finland

DATE OF REVISION OF THE TEXT: 07-FEB-2011

V-2011

PRINT & PACKAGE
JAKKO-TAARA OY
Tel: +358 (0)207 917 600
Fax: +358 (0)207 917 649
Email: repro.ji@jt.fi
www.jt.fi
JAKKO-TAARA OY
RUSKONIE 4
20200 TURKU
FINLAND
JT Job No.: 216967
Item No.: 83261634
Drawn up/OS: Miia Vainio
Date: 24/10/2011
Client: Bayer Oy
Colors: 2/0 321, musta
Pharmacode: 2469
Folding:
Dimensions (mm): 490 x 280 mm
Folded (W x L mm): 70 x 280 mm

Font / Font size / Spacing:
Helvetica Neue / 10 pt

☐ New proof required

☐ Proof approved

Colours are referential.
Please add changes as clearly as possible.
NOTE: We will charge changes that are against manuscript.
JT is not responsible of mistakes in approved proof.

Approved for printing / Date