





Base : 42mm x 130mm

*All text must be 2mm away from diecut outline and glue area.


Please change the following settings in Acrobat Reader to view this file correctly.
For PC: Edit > Preferences > Categories: Page Display > Use Overprint Preview > Select "Always"
For MAC: Adobe Reader > Preferences > Categories: Page Display > Use Overprint Preview > Select "Always"

RE-TYPESET ARTWORK. Due to the nature of the files received, this artwork has been re-typeset.

Varnish Thumbnail 50%



80% Barcode Magnification Truncatch



 inova pharmaceuticals		ARTWORK		PLEASE NOTE This HI RES PDF is not to be used as Final Artwork Supplier information may be present on the packaging EMO/Printing Vendor: Please share the Artwork PDF to Printer for checks on: - Colours - No. of colours, Sequence, CMYK/PMS - Text - Line Thickness, Font Size - Dieline - Size, Eyemark, Print Direction - Coding - Size, Location, VFA, Format Bkg Colour Moving forward when possible any additional information and technical check required by EMO on the Artwork BEFORE APPROVING.	
Job No. : JB-03122 UIC : A8937 SKU : 100500		13.08.24 - V0 - Artwork build - Koragot 11.03.25 - V1 - Amend - Suphapon			
A8937_NUELIN SYRUP 100ML MY Component : FAF LABEL-LEAFLET Country: MY					
Dieline Ref: XXXXXX XXXXXX EMO: ENSIGN Printer: XXXXXX XXXXXXXXXXXX		G247 File Name: 104988_A8937_NUELIN SYRUP 100ML MY_MY_LL_V1		Colour 100%	
		PMS 318 C PMS 285 C		Varnish Guide	