

Technical Validation 23 / 07 / 2020 Version 1		C.C.
Please send open illustrator CC or PDF files for final approval		
Text needs to be converted to outlines		
Use pantone colours solid coated, or CMYK (except from the colour black which must be the one from CMYK)		
Leave 1.5 mm bleed around the knife		
 Remedica FOR A HEALTHIER WORLD		

All the information added by us on the technical validation should be added by the client before sending the artwork for print proof! Once the instructions have been carried out please remember to delete the comments in red colour



CLIENT / PROJECT 	WORK NAME	SIZE	COLOUR		OTHERS	SCALE 1:1
	Ingefitinib_Packaging_Malaysia_Design		 CMYK Printing  Intega Logo (C 67 M 5 Y 9 K 0)  (C 99 M 35 Y 0 K 85)  Ingefitinib Logo  (C 0 M 0 Y 0 K 65)  (C 21 M 0 Y 81 K 0)	 (C 21 M 0 Y 81 K 0)  (C 2 M 0 Y 44 K 0)  (C 0 M 2 Y 0 K 0)   Text  (Grayscale K 100)  (Grayscale K 70)  (C 0 M 100 Y 100 K 0)	 - Die-Cut Line (Do not Print)  - Fold Line (Do not Print)	ARTBOARD A2
						DATE
						ARTWORK BY