




DIMENSION: 39x92 mm (4 tablets per blister)

STANDARD TECHNICAL INFO

- Please do not change the technical info colors - *Cutting*.
- Do not use hidden layers.
- Overprint color usage must be indicated in wording outside the AW.
- **All text on the AW must be vectorized (converted into outlines / curves).**
- All pictures used in AW must be in high resolution and embeded into the file.

DATING OPTIONS:

2-line	- Please do not move the dating position on the AW.
LOT:	- Lek standard dating is LOT: ZZ9999 and MNF/EXP: MM YYYY.
EXP:	

2D DATAMATRIX

- For Datamatrix code creation (combined with dating info), any additional info must be indicated outside the AW (e.g. CIP number).