

FNB4984

ยาเส้นอายุ
2000 mL
(APPROX 80 mL EXCESS)

Vantive

EXTRANEAL Peritoneal Dialysis Solution with 7.5% Icodextrin

EACH 100 mL CONTAINS

7.5 g ICODEXTRIN 538 mg SODIUM CHLORIDE USP 448 mg SODIUM LACTATE
25.7 mg CALCIUM CHLORIDE USP 5.08 mg MAGNESIUM CHLORIDE USP pH 5.2 (5.0 TO 6.0)
mmol/L SODIUM 132 CALCIUM 1.75 MAGNESIUM 0.25 CHLORIDE 96 LACTATE 40

OSMOLARITY 284 mOsmol/L (CALC)

STERILE NONPYROGENIC STORE BELOW 30°C SEE INSERT

FOR INTRAPERITONEAL ADMINISTRATION ONLY NOT FOR INTRAVENOUS USE

CAUTIONS SQUEEZE AND INSPECT INNER BAG WHICH MAINTAINS PRODUCT STERILITY
DISCARD IF LEAKS ARE FOUND DO NOT USE UNLESS SOLUTION IS CLEAR
DISCARD UNUSED PORTION KEEP OUT OF REACH OF CHILDREN

MANUFACTURED BY

VANTIVE MANUFACTURING PTE. LTD.

2 WOODLANDS INDUSTRIAL PARK D STREET 2 SINGAPORE 737778

DIRECTION TO BE USED AS DIRECTED BY PHYSICIAN

FOR HONG KONG ONLY: PRESCRIPTION DRUG 處方藥物

IMPORTED BY

VANTIVE HEALTHCARE (THAILAND) CO., LTD., BANGKOK

FOR INDONESIA ONLY: HARUS DENGAN RESEP DOKTER

IMPORTED & MARKETING BY PT KALBE FARMA TBK, BEKASI – INDONESIA

FOR MALAYSIA ONLY:

CONTROLLED MEDICINE. JAUHI DARIPADA KANAK-KANAK

PRODUCT REGISTRATION HOLDER:

VANTIVE SDN BHD

B-21-3, LEVEL 21, THE ASCENT, PARADIGM, NO. 1, JALAN SS 7/26A, KELANA JAYA,
47301 PETALING JAYA, SELANGOR

ULTRABAG CONTAINER PL-146
SIN11220P MAL06021264AZ
HK-46627 BRU15011744NP
INDONESIA DKI 0915600242A1
THAI REG. NO. 2C 16/51 (N)

請儲存於攝氏 30 度以下
用法：請依照醫生指示使用

"萬益特" 愛多尼爾腹膜透析液

衛署藥輸字第 023687 號
請詳閱說明書

本藥限由醫師使用
萬益特醫療用品股份有限公司

批號及保存期限詳如包裝所示
台北市信義區忠孝東路五段
68 號 20 樓 A1 室

Vantive, EXTRANEAL, UL-TRABAG & PL-146 are trademarks of Vantive Health LLC or its affiliates.