

Document Number: 15000

KONTRAK KKM

45mm

Oxacil Granules 125mg/5ml **60ml**

Each 5ml contains: Cloxacillin Sodium eq. to Cloxacillin 125mg, Sodium Benzoate 0.2%w/v

Please refer package insert for full information.

Route of administration : Oral

CONTROLLED MEDICINE Not to be used in patients with known hypersensitivity to penicillin.

Keep container tightly closed. Store in a dry place below 30°C. Protect from light. Store in refrigerator after reconstitution. Discard 7 days after reconstitution.

Keep out of reach of children. Jauhi daripada kanak-kanak.

Instruction for reconstitution: Add water to the mark on bottle (60ml). Shake well and make up to the mark if needed. Share well before taking each dose.

Batch No. :
Mfg. Date :
Exp. Date :
Reg. No. : MAL19860043AZ
15000 KONTRAK KKM

Product Registration Holder and Manufacturer:
DUOPHARMA (M) SDN. BHD.
Lot. 2599, Jalan Seruling 59, Kawasan 3, Taman Klang Jaya,
41200 Klang, Selangor Darul Ehsan, Malaysia.

76mm



Black



Pantone 348 C



Pantone Red 032 C



Pantone 375 C (Transparency 40%)

As per printer advised 29.3.2016



0c 100m 100y 0k
PANTONE 2035 C



0c 100m 100y 50k
PANTONE 7624 C

Scale 150%

45mm

Oxacil Granules 125mg/5ml **60ml**

Each 5ml contains: Cloxacillin Sodium eq. to Cloxacillin 125mg, Sodium Benzoate 0.2%w/v

Please refer package insert for full information.

Route of administration : Oral

CONTROLLED MEDICINE Not to be used in patients with known hypersensitivity to penicillin.

Keep container tightly closed. Store in a dry place below 30°C. Protect from light. Store in refrigerator after reconstitution. Discard 7 days after reconstitution.

Keep out of reach of children. Jauhi daripada kanak-kanak.

Instruction for reconstitution: Add water to the mark on bottle (60ml). Shake well and make up to the mark if needed. Share well before taking each dose.

Batch No. :
Mfg. Date :
Exp. Date :
Reg. No. : MAL19860043AZ
15000 KONTRAK KKM

Product Registration Holder and Manufacturer:
DUOPHARMA (M) SDN. BHD.
Lot. 2599, Jalan Seruling 59, Kawasan 3, Taman Klang Jaya,
41200 Klang, Selangor Darul Ehsan, Malaysia.

76mm