

200% Enlarged

Kauliv 20 µg / 80 µL

solution for injection
teriparatide

subcutaneous use

Store in refrigerator at 2°C to 8°C

Keep out of reach of children

Jauhkan ubat daripada kanak-kanak

Use only with Kauliv Pen

3 mL

LOT:

EXP:

NO VARNISH ZONE

WINDOW SIZE: 40 x 14 mm

1200002693

UNWIND DIRECTION OF LABEL ROLL

Kauliv 20 µg / 80 µL

solution for injection
teriparatide

subcutaneous use

Store in refrigerator at 2°C to 8°C

Keep out of reach of children

Jauhkan ubat daripada kanak-kanak

Use only with Kauliv Pen

3 mL

LOT:

EXP:

NO VARNISH ZONE

WINDOW SIZE: 40 x 14 mm

1200002693

50 x 40 mm
ACTUAL SIZE