










Artwork Code : XXXXXXXXXX-XX				Artwork Req. No: -					
Quality of Paper/Board : Chromo				Size of Artwork (In mm) : 53x53					
Quality of Gum : 19 GSM				GSM of Paper/Board : 80 GSM					
Colour Code :	BLACK	P 575	Pantone	Pantone	Pantone	Pantone	Pantone	Pantone	
									
Barcode Information:									
Barcode Scan Report:									
Packing: 200 mL PLB (Euro head Closed)				Plant Location : MF 02					
Country: Malaysia				Language : English, Malay					
Ref. Code Creation/Blockage Note:				Artwork Control Key No.:					
"Controlled Copy" Holder (1)			(2)		(3)		(4)		(5)


Otsuka
100 ml

**Sterile, non-pyrogenic, single dose container:
Sodium chloride solution for intravenous infusion BP
(0.9 % w/v)**

NS

Caution : Even invisible damage to container caused during transit, storage or handling may result in contamination. Do not use, if container is found leaking upon squeezing or solution is not clear. Solution containing visible solid particles must not be used.

Each 100 ml contains :
Sodium Chloride BP 0.9 g
Water for Injections BP q.s.

M. L.: G/28D/LV/P5

B. No.:



Mfg. Dt.:

Exp. Dt.:

MAL No.:

Discard all unused contents.
Dosage : As directed by the Physician.
Store at a temperature not exceeding 30°C.
To be used with non-pyrogenic I.V. administration set with aseptic technique.
For IV infusion use only.
JALUHI UBAT DARIKANAK-KANAK
KEEP MEDICINE OUT OF REACH OF CHILDREN
CONTROLLED MEDICINE/UBAT TERKAWAL

Otsuka
Manufactured by:
Otsuka Pharmaceutical India Private Limited
Survey No. 159 to 201 & 208 to 210,
Village-Vasana-Chacharwadi, Tal - Saranad,
Dist. Ahmedabad - 382 213, India.


Product Registration Holder, Importer & Distributor:
BioCare Pharmaceutical (M) Sdn. Bhd.
No.285, 28-1 & 28-2,
Pusat Dagangan P/2C,
Jalan P/5 S/28 46150 Petaling Jaya,
Selangor, Malaysia.

XXXXXXXXXX

15.00 mm
Batch Printing Area

2.00 mm
Batch Printing Area

15.00 mm
Batch Printing Area



"Artwork number will be generated at the time of commercial process"