



Budesonide Inhaler

200 mcg / dose

Budecort 200

300
METERED
DOSES

Cipla



Budecort 200

Each actuation delivers
Budesonide BP 200 mcg
Suspended in
Propellant HFA-134a q.s.

Directions for use:
Read instruction leaflet enclosed before use. Use only as directed by the physician.

Route of administration:
Oral Inhalation

It is dangerous to exceed the recommended dose.

Warning: Pressurised canister, keep away from sunlight and heat. Do not puncture, break or burn even when apparently empty.



Budesonide Inhaler

200 mcg / dose

Budecort 200

300
METERED
DOSES

Cipla



Budecort 200

Store below 30°C.
Do not freeze.
Keep medicines out of reach of children / Jauhi ubat dari kanak-kanak. Please refer enclosed leaflet for more information.

Controlled Medicine/ Ubat Terkawal
Reg. No: MAL08051745ACZ
M.L. 341/L

Registration Holder in Malaysia
CIPLA MALAYSIA SDN BHD (Company No. 1039096-P)
Suite 1101, Amcorp Tower, Amcorp Trade Centre, 18 Persiaran Barat, 46 050 Petaling Jaya, Selangor, Malaysia.

Mfd. for **CIPLA LTD.**
At Medispray Laboratories Pvt Ltd. 344/345 Kundaim Indl. Estate Kundaim, Goa-403115, India.

B. No.
MFD.
EXP.

xxxxxxxx

← Unvarnished Area

PACKAGING DEVELOPMENT

Product Name: Budecort 200 Inhaler		Material No.: xxxxxxxx	Version: 01	Item: Carton	Co-ordinator: Shruti	Artist: Vasundhara	Date: 15-02-2018																																
Colours: BLUE WOOL TEST VALUE 5-8 (LIGHT FASTENING DATA) ■ Pantone 226 C ■ Pantone 620 C ■ Pantone 356 C ■ Pantone 485 C				Ink: Oil Base Colours NBCL Series from SIEGWERK OR Oil based inks from DIC. Varnish: Aqua-Water based OP Ink-High Gloss High Slip(10-603787-2,2660) from SIEGWERK OR HIV-1611-High Gloss Aqua Duct OPV from MICRO INKS																																			
Spectro-Densitometer Delta-E reading (ΔE) for colour: NOT MORE THAN dE2.5		Fonts: ----			Software: Illustrator CC																																		
Design: Cartonator (CAM)				Links: ----																																			
Actual Size: 50 x 35 x 96 mm		Pharmacode: --_STD		Grain Direction : Perpendicular to Crease		Supersedes/Reference: 21049634																																	
Tuck flap: 15 mm		Side / Collar flap: 24 mm		Overlap: No		Caliper (Thickness) for Board: NLT 430 micron																																	
Material: 300 GSM CFB Board (Cyber XL from ITC /Tuffcote from J.K. PAPER) Overprinting area to be left unvarnished.				Varnish: Scuff Proof Aqua		Artwork Print Size: <input type="checkbox"/> actual <input type="checkbox"/> scaled																																	
Path: PC:/D:/Vasundhara/Malaysia/Budecort/xxxxxxxx Budecort 200 Inhaler 300MD Carton Malaysia.ai																																							
<ul style="list-style-type: none"> Instructions/Remark: ----- Any deviation must be brought to the notice of packaging development co-ordinator immediately. For any clarification, please contact packaging development co-ordinator immediately. NO CHANGES IN ARTWORK SHOULD BE DONE BY THE PRINTER The printer should verify the e-proof against the approved artwork before submitting for approval and the e-proof should have printer details. 				<table border="1"> <tr><th>Checked by</th><th>Artist</th><th>Cordinator</th><th>file loaded in Server</th><th>Section Head</th></tr> <tr><td>Pharma Code</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> <tr><td>2D Code</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> <tr><td>QR Code</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> <tr><td>Bar Code</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> <tr><td>Artwork</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> <tr><td>Spell check</td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td></td><td></td></tr> </table>	Checked by	Artist	Cordinator	file loaded in Server	Section Head	Pharma Code	<input type="checkbox"/>	<input type="checkbox"/>			2D Code	<input type="checkbox"/>	<input type="checkbox"/>			QR Code	<input type="checkbox"/>	<input type="checkbox"/>			Bar Code	<input type="checkbox"/>	<input type="checkbox"/>			Artwork	<input type="checkbox"/>	<input type="checkbox"/>			Spell check	<input type="checkbox"/>	<input type="checkbox"/>		
Checked by	Artist	Cordinator	file loaded in Server	Section Head																																			
Pharma Code	<input type="checkbox"/>	<input type="checkbox"/>																																					
2D Code	<input type="checkbox"/>	<input type="checkbox"/>																																					
QR Code	<input type="checkbox"/>	<input type="checkbox"/>																																					
Bar Code	<input type="checkbox"/>	<input type="checkbox"/>																																					
Artwork	<input type="checkbox"/>	<input type="checkbox"/>																																					
Spell check	<input type="checkbox"/>	<input type="checkbox"/>																																					