




PACKAGING DEVELOPMENT

Product Name: Ipracip Inhaler 200MD		Material No.: XXXXXX	Version: 01	Item: Label	Co-ordinator: Vinod	Artist: Mukesh	Date: 30-03-2020	
Colours: BLUE WOOL TEST VALUE 5-8 (LIGHT FASTENING DATA) Pantone 186 C Pantone 356 C Pantone Orange 021 C Black					Ink: Aqua (UNI Paper series) from SIEGWERK Varnish: Aqua Water based OP Ink - High Gloss High Slip (10-603787-2.2660) from SIEGWERK			
Spectro-Densitometer Delta-E reading (ΔE) for colour: NOT MORE THAN dE2.5		Design: Roll Form			Reference: 217259		Software: Illustrator CC	
Fonts: ----				Links: -				
Actual Size: 68 x 35 mm	2 D Code:-	Grain Direction : Parallel to length			Unwinding Direction: As indicated in Artwork	Screen : # __		
Material: 80 gsm Fasson FasPrint NG paper/Adhesive Fasson HITAC/ PET Release liner Fasson PET23 (PTA2450 from Avery Dennison) with adhesive NLT 15 gsm and PET release liner NLT 28 gsm or 23 micron ±10 %					Varnish: Spot Scuff Proof Aqua		Artwork Print Size: ☐ actual ☐ scaled	
Path: PC D : Mukesh\Vinod\Malaysia\Ipracip : Ipracip Inhaler 200MD Label Malaysia.ai								
• Instructions / Remark: ---- • Any deviation must be brought to the notice of packaging development co-ordinator immediately. • For any clarification, please contact packaging development co-ordinator immediately. NO CHANGES IN ARTWORK SHOULD BE DONE BY THE PRINTER. • The final sample set (shade card, multi-up sheet and samples as per CIPLA Standards) should be sent FOR CIPLAS Approval and it should be signed and dated, and also be accompanied with one multi up which is approved by Supplier's QA/QC with sign and date and should have Stamped with "Copy matter and colour checked for all the multi ups" on the multi up sheet.				Checked by	Artist	Cordinator	file loaded in Server	Section Head
				Pharma Code	☐	☐		
				2D Code	☐	☐		
				QR Code	☐	☐		
				Bar Code	☐	☐		
				Artwork	☐	☐		
Spell check	☐	☐						