

Figure 1 illustrates the dimensions of the printed matter area for the Rosuvastatin Film-Coated Tablets USP 20 mg packaging. The figure is divided into two main sections: a top view and a side view.

Top View Dimensions:

- Width: 43.965 mm
- Depth: 43.00 mm

Side View Dimensions:

- Height: 224.00 mm
- Width: 87.7 mm - PRINTED MATTER AREA

The packaging is labeled with the following information:

- Top View:**
 - UMEDICA
 - Mfg. Lic. No. : G/25/491
 - EXP MM/YYYY
- Side View:**
 - UMEDICA
 - Mfg. Lic. No. : G/25/491
 - EXP MM/YYYY

The image displays a vertical strip of five identical white rectangular labels, each representing a box of SYNERRV ROSUVASTATIN 20 Rosuvastatin Film-Coated Tablets USP 20 mg. The labels are arranged in a column, separated by thin white gaps. Each label has a clean, minimalist design with a white background. At the top of each label is a small, stylized globe icon. Below the icon, the text 'SYNERRV ROSUVASTATIN 20' is printed in a bold, black, sans-serif font. Underneath this, 'Rosuvastatin Film-Coated Tablets USP 20 mg' is written in a slightly smaller, black, sans-serif font. At the bottom of each label, the 'UMEDICA' logo is displayed, consisting of a globe icon followed by the word 'UMEDICA' in a bold, black, sans-serif font. To the right of the logo, the text 'Mfg. Lic. No. : G/25/491' is printed in a smaller, black, sans-serif font. The entire strip is set against a dark, textured background.

Format																					
Format Name		SPECIFICATION FOR ARTWORK APPROVAL																			
Format Number		FM/PDD/002-01-07				SOP Number				PDD/002											
New Artwork Code No. : 205570		Old Artwork Code No.: NA				Artwork Effective Date:															
Product Name		: SYNERRV ROSUVASTATIN 20 (Rosuvastatin Film-Coated Tablets USP 20 mg)				Company Name/Country : Synerrv / Malaysia															
Component Type		: Blister Foil				Pack Style : 10T - Alu Alu Blister Foil															
Description		: New Development Blister Foil Artwork as per Export-Order/Malaysia Country requirement, as per Change Control No. CC/PDD24/056.																			
Specifications : LxWxH Layout/Drawing No.: (PET122399E675-R1) Site : Umedica Artwork Language: English																					
Full Size / Size		: NA				Type+ Blister +Strip Size+ Adhesive+ Tube Type				: 91x38x4.2 (TB) MM											
Size –Catch Cover/ Foil +BOPP		: 224 MM ± 1 MM (223 to 225 MM)				Sachet Size/Pocket Size/Fold Size/Release Paper :				NA											
Space between Labels		: NA				Transparent Paper Thickness				: NA											
GSM of Paper+ GSM of Aluminium Foil :		Between 63.5 & 71.6 g/m²				Type of Paper / GSM				: NA											
Thickness/Thickness of Foil		: 0.025 MM (24µ to 27µ)				Type of Foil+ Fold+ Coating GSM				: Between 6.0 & 8.0 g/m²											
Cylinder		: 1 Nos				Lamination / Varnish / Material Grade				: NA											
Repeat/Tolerance : 43.965 MM ± 0.025 MM						Partition / Plate Size / Total GSM /				: NA											
(20 Repeats = 878.8 ~ 879.8 MM) [Eyemark]						Sealing / Embossing / Core Diameter				: NA											
GSM/BF of Printed Shipper		: NA				Paper Grade / Flute Type / PLY of Shipper				: NA											
(Top/Middle/Bottom and Flute layer		: NA																			
Printing Colour		: P-485 C [Red Shade]								or as per Sample											
Approved by Export/Marketing : [Approval attached with ADF]																					
CHECKLIST FOR CHECKING ARTWORK & SPECIFICATION																		Artwork Approved by Sign with Date			
Brand Name and Generic (P/P/USP/In-house copy)	Dosage Form and strength	Prefix Rn/ XRx	Mfg. Lic. No. /Neutral Code	Regulatory Approval	Packing (Qty) / Presentation	Dosage details if required	Size / Design	Composition label claim/ Formulation	H, H+, C+ or X, ac for local product / Fed line for genl products	Item Code	QR Code/ Barcode (1D 2D)/ Plasma code requirement	Material and printing specifications	Storage Condition as per Buyer Registration	Registration No. / Other requirement	Mfg. Name and Address (Rel. to the insert)	Manufactured by/ Distributors Name	Route of Admin. (in Multiple) for Overprinting	Adequate Space (Bath or embossing) for Overprinting	Prepared by PDD Department	Approved by PDD Department	
																				QC Department	
																				RA Department	
																				PPIC Department	Production Packing Department
																				QA Department	
Note to Printer:																					
• After receive the CDR File, Submit the print proof for approval.																					
• The overprinting area should be unvarnished.																					
• Machine proof must be verified against the approved hard copy in Light, Standard and Dark shade should be certified and signed by an authorized person.																					
• Any deviation must be brought to the notice of the Packaging Development Department immediately in the writing.																					
• For any clarification, please contact Packaging Development Department.																					