


Item Number: 0755504	Version: 1	Swatches	
Profile: 0402054_4_UC_Drop	Technical Specification: 0402054	Visible	G M Y K
Artwork Dimensions: Size: 145 x 30 x 200 mm			
Graphics House: SGK Richmond		Visible / Not Printed*	
Date: 01SEP2021		Dieline	Technical Info
OPTIONAL: Artwork Approver outside the Shire Management System:			
Role:	Name:	Signature:	Date:
			
*Color not present in printed material		Body Text Size: 7.0 pt	Code 1TF 2/5 755504
		Datamatrix Code NA	RSS / GS1-128 Code NA
		Pharmacoce NA	Code 3 of 9 NA
			GTIN / EAN-13 Code 0365412649223

4 mL



Frozen solution



0 853412 649223

Reorder number: 5500172

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe



4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL



Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL



Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe

4 mL

Frozen solution

4 mL

Frozen solution

TISSEEL 4 mL

[Fibrin Sealant]

Frozen solution
with Pre-filled PRIMA Syringe