

Carton: 2010-22-006-28-316
LIFT - Current 4

Lacquerform 30883-1
Lacquer free area: ■

Braille No.: XXXX-XX-X
Placement: Front

Colour:
PMS 626C + PMS 134C +
PMS 280C + PMS 555C + Black

Manuf./Expiry/Batch:

Use a new needle for every injection

1x1 mL

once weekly

700 units/mL

Awqli®
FlexTouch®

insulin icodec
solution for injection in pre-filled pen

22-006-28

Keep out of the sight and reach of children/luar ni ubat dari kanak-kanak
Read the package leaflet before use
Use only clear, colourless solution
Before first use: Store in a refrigerator (2°C - 8°C)
Do not freeze
After first opening: Store below 30 °C or in a refrigerator (2°C-8°C)
Use within 12 weeks
Keep the cap on the pen in order to protect it from light
Discard the needle safely after each injection
Single patient use only
1 mL solution contains 700 units of insulin icodec (equivalent to 26.8 mg). Each pre-filled pen contains 700 units of insulin icodec in 1 mL solution
Excipients: Glycerol, metacresol (1.00 mg/mL), phenol (5.65 mg/mL), zinc acetate, sodium chloride, hydrochloric acid (for pH adjustment), sodium hydroxide (for pH adjustment) and water for injections. See leaflet for further information
Controlled Medicine/Ubata Terkawal



SV
perd
:SV
M
ped

once weekly

insulin i codec
solution for injection in pre-filled pen

700 units/mL

Awikli®
FlexTouch® 700 units/mL

solution for injection in pre-filled pen
insulin icodec
subcutaneous use

once weekly MALXXXXXXXXXAZ

1 x 1 mL pre-filled pen (700 units) with 9 disposable needles

Awigili®
FlexTouch®

 The pen shows the dose
One step equals 10 units



Open here

8-9529-07-301-1.indd 1

3/03/2025 16:19:32



Project: PM2025-00226
Document: 8-9529-07-301-1_v11.pdf
Document description:
Document Lifecycle State: Mock-up
Version: 1
Version ID: 00002_0000936116
Downloaded: Mar 4, 2025 at 10:40 AM GMT+08:00
By: Yunny, Tung (TPYG@NOVONORDISK.COM)
Approval state: **Draft**

This is a representation of an electronic record that was signed electronically and this page is the manifestation of the electronic signature.
