

Code centre line  
Code: 100% Direction →  
Length: Max 29 mm (100%)

Regulatory Operations  
Carton:  
2010-22-503-57-302  
Current 1.0  
Lacquer free area: ■  
Lacquerform: 30100-1  
Braille No.:  
Placement: Front  
Colour: PMS 280C  
+ PMS 289C  
+ PMS 471C + Black

2 —  
8 —



1 ml suspension contains 100 IU (3.5 mg) of insulin human (rDNA) (30% as soluble insulin and 70% as isophane insulin), zinc chloride, glycerol, metacresol 1.5 mg, phenol 0.65 mg, disodium phosphate dihydrate, sodium hydroxide, hydrochloric acid, protamine sulphate and water for injections

**Keep out of the reach and sight of children**  
**Read the package leaflet before use**  
**Resuspend according to instructions**  
**Single patient use only**  
**Store in a refrigerator (2°C - 8°C)**  
**Do not freeze**  
**Keep the pen cap on in order to protect from light**  
**After first opening: Use within 6 weeks when stored below 30°C. Do not refrigerate**  
**Discard the needle after each injection**

Jauhi ubat dari kanak-kanak  
This is a STERILE product  
Controlled Medicine

**Product Registration Holder:**  
Novo Nordisk Pharma (M) Sdn Bhd  
Menara 1 Sentrum, Level 16  
No. 201, Jalan Tun Sambanthan  
50470 Kuala Lumpur

**Manufactured by:**  
Novo Nordisk A/S,  
Novo Allé  
DK-2880 Bagsværd  
Denmark

© 2016  
Novo Nordisk A/S



22-503-57

8-9653-00-310-1

**Mixtard® 30**  
FlexPen®  
100 IU/ml

100 IU/ml  
FlexPen®  
30

100 IU/ml  
FlexPen®  
30

5x3 ml  
**Mixtard® 30** 100 IU/ml  
100 IU/ml  
Suspension for injection  
Pre-filled pen, FlexPen®  
Insulin human (rDNA)

Subcutaneous use  
Designed to be used with NovoFine® or NovoTwist®  
disposable needles up to a length of 8 mm  
Needles are not included  
MAL13115134ASZ



**Mixtard® 30**  
FlexPen®  
100 IU/ml  
Manuf./Expiry/Batch:  
[Redacted]


Regulatory Operations

Carton: 2010-22-476-07-301  
Current 4.0

Lacquerform 30123-1  
Lacquer free area: ■

Braille No.: XXXX-XX-X  
Placement: (Front)

Colour: PMS 280C +  
PMS 471C



2,9 mm  
8 mm  
1 mm

Code centre line  
Code: 100% Direction →  
Length: Max. 29 mm (100 %)

2 mm  
5 mm  
2 mm