

INFUSOL® HMK
DARROW'S SOLUTION
INTRAVENOUS INFUSION BP

Sterile and Non-Pyrogenic

The solution contains:

Sodium Chloride BP	4.00 g/L
Potassium Chloride BP	2.60 g/L
Sodium Lactate USP	5.94 g/L
Na ⁺	121.00 mmol/L
K ⁺	35.00 mmol/L
Cl ⁻	103.00 mmol/L
HCO ₃ ⁻ (as lactate)	53.00 mmol/L
Osmolarity:	312.00 mOsm/L

FOR INTRAVENOUS USE

CAUTION : Single dose container. Discard all unused contents. Do not use if leakage is detected. If any visible solid particles, growth or turbidity appears during storage, the product should not be used.

STORAGE : Do not store above 30°C

*Keep out of reach of children
*Jauhkan daripada kanak-kanak

500mL

MAL20014358AZ

Controlled Medicine / Ubat Terkawal

Manufacturer / Product Registration Holder:
AIN MEDICARE SDN. BHD.
 Jalan 6/44, Kaw. Perindustrian Pengkalen Chepa 2
 16100 Kota Bharu, Kelantan, MALAYSIA.

MS 2424:2012
1 062-09/2013

DARROW'S SOLUTION

Remarks: All the following information will be provisioned on the label;

- 1) Hologram sticker
- 2) QR/Bar code
- 3) Registration number of other countries.

APPROX

APPROX

100

100

200

200

300

300

400

400

INFUSOL® HMK

**DARROW'S
SOLUTION**

INTRAVENOUS INFUSION BP

Sterile and Non-Pyrogenic

The Solution Contains:

Sodium Chloride BP	4.00 g/L
Potassium Chloride BP	2.60 g/L
Sodium Lactate USP	5.94 g/L
Na ⁺	121.00 mmol/L
K ⁺	35.00 mmol/L
Cl ⁻	103.00 mmol/L
HCO ₃ ⁻ (as lactate)	53.00 mmol/L
Osmolarity:	312.00 mOsm/L

FOR INTRAVENOUS USE

CAUTION: Not for use in the treatment of lactic acidosis.
Discard any portion not used when the seal is first broken. Before use check for minute leaks by squeezing bag firmly. Discard if leakage is detected. If any solid particles, growth or turbidity appears during storage, the product should not be used.

STORAGE: Do not store above 30°C

* Keep out of reach of children.
* Jauhkan daripada kanak-kanak.

Controlled Medicine / Ubat Terkawal

500 mL MAL20014358AZ



Manufacturer/Product Registration Holder:
AIN MEDICARE SDN. BHD.
Jalan 6/44, Kaw. Perindustrian Pengkalan Chepa 2,
16100 Kota Bharu, Kelantan, MALAYSIA.



MS 2424:2012
1 082-09/2013

Remarks :

The following information will be provisioned on the label;

1)Registration number of other countries

Manufacturing date, Expiry date and Batch No. will be printed during production.

APPROX.
100
200
300
400
500
600
700
800
900

INFUSOL® HMK

DARROW'S

SOLUTION

INTRAVENOUS INFUSION BP

Sterile and Non-Pyrogenic

The solution Contains:

Sodium Chloride BP	4.00 g/L
Potassium Chloride BP	2.60 g/L
Sodium Lactate USP	5.94 g/L
Na ⁺	121.00 mmol/L
K ⁺	35.00 mmol/L
Cl ⁻	103.00 mmol/L
HCO ₃ ⁻ (as lactate)	53.00 mmol/L
Osmolarity:	312.00 mOsm/L

FOR INTRAVENOUS USE

CAUTION: Not for use in the treatment of lactic acidosis.
Discard any portion not used when the seal is first broken. Before use check for minute leaks by squeezing bag firmly. Discard if leakage is detected. If any solid particles, growth or turbidity appears during storage, the product should not be used.

STORAGE: Do not store above 30°C

*Keep out of reach of children
*Jauhkan daripada kanak-kanak

Controlled Medicine / Ubat Terkawal

1000 mL **MAL20014358AZ**



Manufacturer/Product Registration Holder:
AIN MEDICARE SDN. BHD.
Jalan 6/44, Kaw. Perindustrian Pengkalan Chepa 2,
16100 Kota Bharu, Kelantan, MALAYSIA.



APPROX.
100
200
300
400
500
600
700
800
900

Remarks :

The following information will be provisioned on the label;

1)Registration number of other countries

Manufacturing date, Expiry date and Batch No. will be printed during production.