




Black

Outer & Dotted Lines Not to be Printed.

Blister size 95 x 42 mm

Foil Width 101 mm


Continuous Printing

Layout No: BFBQS-B5929/2

UNWINDING DIRECTION

Foil unwinding direction on Machine

		Product Name	Component	Item Code	Date & Time
		BISODAC-HF 2.5 mg	Foil (10's) Push Through Blister	P0116139	12.12.2018 & 4:50 pm
Customer / Country		Version No.	Reason Of Issue	Reviewed / Approved by	
Malaysia		00	Commercial		
Team Leader	Lova Raju	Dimensions	No. of Colours : 01		
Initiator	Bharath	95 x 42 mm (101 mm)			
Artist:			Pharmacode	Sign / Date	
		NA			
Additional Information :		WOS-18-000919			



■ Black

Blister size 103 x 58 mm

Total Width 109 mm

L No: BFBQS-A3516/4

→ UNWINDING DIRECTION → Foil unwinding direction on Machine

Outlines not to be printed
Blisters are with perforation

Dotted lines are not to be printed

 AUROBINDO Packaging Development		Product Name	Component	Item Code	Date & Time
		BISODAC-HF 2.5	Foil 10's	P0116137	12.12.2018 & 10.00 am
		Country/Customer	Version No.	Reason of Issue	Reviewed / Approved by
		Malaysia	00	Submission	
Team Leader	Lovaraju	Dimensions (mm)	Colours		
Initiator	Bharath	103 x 58 mm (109 mm)	01		
Artist	Sree Designers				
Additional Information :					
				Sign / Date	