

SUSPENSION

CEPHALEXIN-250



CEPHALEXIN-250

SUSPENSION

Each 5ml contains:
Cephalexin Monohydrate
equivalent to
Cephalexin 250mg

60ml

HOW TO PREPARE YOUR MEDICINE



TO COMPLETE THE COURSE OF MEDICINE

Barcode



CEPHALEXIN-250

SUSPENSION

Each 5ml contains:
Cephalexin Monohydrate
equivalent to
Cephalexin 250mg

60ml

60ml when reconstituted.
For oral administration only.

PRECAUTION:
Should be used with caution in patients with penicillin sensitivity.

Keep container well closed. Store below 30°C until reconstituted. After reconstitution, refrigerate between 2°C to 8°C and use within 7 days.

Please refer to enclosed package insert for more information.

**KEEP OUT OF REACH OF CHILDREN
JAUHI DARI KANAK-KANAK**

Controlled Medicine | SHAKE WELL BEFORE USE
GONCANG SEBELUM DIGUNAKAN

MAL16035027AZ

Manufacturer & Product Registration Holder :



Kotra Pharma (M) Sdn. Bhd.
No.1, 2 & 3, Jalan TTC 12, (50082-V)
Cheng Industrial Estate,
75250 Melaka, Malaysia.

50x62x120

BATCH NO.
MFG DATE :
EXP DATE :

SUSPENSION

CEPHALEXIN-250



CEPHALEXIN-250

SUSPENSION

Each 5ml contains:

Cephalexin Monohydrate

equivalent to

Cephalexin 250mg

100ml when reconstituted.

For oral administration only.

PRECAUTION:

Should be used with caution in patients with penicillin sensitivity.

Keep container well closed. Store below 30°C until reconstituted. After reconstitution, refrigerate between 2°C to 8°C and use within 7 days.

Please refer to enclosed package insert for more information.

**KEEP OUT OF REACH OF CHILDREN
JAUHI DARI KANAK-KANAK**

Controlled Medicine SHAKE WELL BEFORE USE
GONCANG SEBELUM DIGUNAKAN
MAL16035027AZ

Manufacturer & Product Registration Holder :
Kotra Pharma (M) Sdn. Bhd.
No. 1, 2 & 3, Jalan TTC 12, (90082-V)
Cheng Industrial Estate,
KOTRA PHARMA 75250 Melaka, Malaysia.

100ml



CEPHALEXIN-250

SUSPENSION

Each 5ml contains:

Cephalexin Monohydrate

equivalent to

Cephalexin 250mg

100ml

HOW TO PREPARE YOUR MEDICINE

1. SHAKE BOTTLE WELL TO LOOSEN GRANULES
2. FIRST, REMOVE CAP AND DOSE-MASTER™ ADAPTER. THEN, ADD SOME WATER TO MIX WITH THE GRANULES. AFTER THAT, SHAKE AND TOP UP TO 100ML MARK
3. PLACE THE DOSE-MASTER™ ADAPTER AND CAP
4. SHAKE THE BOTTLE WELL TO MIX THE MEDICINE PROPERLY
5. PLACE THE DOSE-MASTER™
6. WITHDRAW THE MEDICINE WITH DOSE-MASTER™ FOR ACCURATE DOSE

TO COMPLETE THE COURSE OF MEDICINE



54 x 65 x128

BATCH NO.
MFG DATE :
EXP DATE :