

Immediate Label For Plastic Bottle of 60 ml:

This Design Also Applicable For Plastic Bottles of 100 ml

<p>Direction for mixing: Shake to loosen powder. Add approximately 30ml freshly boiled and cooled water and shake gently. Add more water until the 60ml mark is reached and shake well.</p> <p>Keep refrigerated (2-8°C) after reconstitution. To be used within 7 days of mixing.</p> <p>FOR ORAL USE ONLY</p> <p>Complete the prescribed course.</p> <p>Not to be used in patients with known hypersensitivity to penicillin.</p> <p>Manufacturer / Product Registration Holder: DYNAPHARM (M) SDN. BHD. 198001011897 (65683-V) 2497, MK.1, Lorong Perusahaan Baru 5, Kawasan Perusahaan Perai 3, 13600 Perai, Pulau Pinang, Malaysia.</p>	<p>Dyna Amoxicillin 250mg/5ml</p> <p>Each 5 ml of the reconstituted suspension contains: Amoxicillin Trihydrate equivalent to Amoxicillin 250 mg</p> <p>Preservative: Sodium Benzoate 0.1% w/v</p> <p>60 ml MAL19985982AZ CONTROLLED MEDICINE</p>	<p></p> <p>SHAKE WELL BEFORE USE</p> <p>JAUHKAN UBAT DARIPADA KANAK-KANAK KEEP OUT OF REACH OF CHILDREN Keep Container Tightly Closed Do Not Store Above 30°C Protect From Light</p> <p>Please refer to the package insert for Indication, Dosage and other information.</p> <p>Batch No.: Mfg. Date: Exp. Date:</p> <p> 9 557490 701050</p> <p>LB1PL A005-10-L25</p>
--	---	--