

The Aluminium Bag (Preprinted With 60 ml x 10 bottles) For Plastic Bottle of 60 ml x 10 Bottles:



[Stick label here]

60 ml x 10 bottles

The Aluminium Bag Is Labelled and/or Stick-on With Product Label (As Proposed To Section D1):

Direction for mixing:
Shake to loosen powder. Add approximately 30ml freshly boiled and cooled water and **shake gently**. Add more water until the 60ml mark is reached and **shake well**.

Keep refrigerated (2-8°C) after reconstitution.
To be used within 7 days of mixing.

FOR ORAL USE ONLY

Complete the prescribed course.

Not to be used in patients with known **hypersensitivity** to penicillin.

Dyna Amoxicillin 250mg/5ml

Each 5 ml of the reconstituted suspension contains:
Amoxicillin Trihydrate equivalent to **Amoxicillin 250 mg**

Preservative:
Sodium Benzoate 0.1% w/v

SHAKE WELL BEFORE USE
JAUHKAN UBAT DARIPADA KANAK-KANAK
KEEP OUT OF REACH OF CHILDREN
Keep Container Tightly Closed
Do Not Store Above 30°C
Protect From Light

Please refer to the package insert for **Indication, Dosage** and other information.

Manufacturer / Product Registration Holder:
DYNAPHARM (M) SDN. BHD.
19801011897 (65685-V)
2497, MK.1, Lorong Perusahaan Baru 5,
Kawasan Perusahaan Perai 3,
13600 Perai, Pulau Pinang, Malaysia.

60 ml
MAL19985982AZ
CONTROLLED MEDICINE

9 557490 701050

LB1PL A005-10-L25