



BCG製品箱マレーシア政府 (5985)



1

3585-04
2019.08



Lot No.:
Mfg. Date:
Expiry Date:

Product Registration Holder:
Eisai (Malaysia) Sdn. Bhd. (18039-D),
Unit 701D, Level 7, Tower D, Uptown 5,
No. 5, Jalan SS21/39, Damansara Uptown,
47400 Petaling Jaya,
Selangor, Malaysia.

-Keep medicine out of reach of children/
Jauhkan ubat daripada kanak-kanak
-Controlled Medicine/ Ubat terkawal

To be used immediately after reconstitution and to discard any unused portion after 6 hours.

FREEZE-DRIED GLUTAMATE
BCG VACCINE (JAPAN)
FOR INTRADERMAL USE
0.5mg BCG × 100 amps.
Storage: +2°C ~ +8°C
MAL19872204AZ
KEMENTERIAN KESIHATAN MALAYSIA
Kod MUS No.07.4201.06

JAPAN BCG LABORATORY
TOKYO, JAPAN



BCG

FREEZE-DRIED GLUTAMATE
BCG VACCINE (JAPAN)
FOR INTRADERMAL USE

0.5mg BCG × 100 amps.
Storage: +2°C ~ +8°C

1 ampoule with 0.5mg BCG (moist weight)
Reconstitute with 1.0ml saline diluent
Do not expose to direct sunlight

JAPAN BCG LABORATORY
TOKYO, JAPAN



JAPAN BCG LABORATORY
TOKYO, JAPAN



JAPAN BCG LABORATORY
TOKYO, JAPAN



JAPAN BCG LABORATORY
TOKYO, JAPAN



FREEZE-DRIED GLUTAMATE BCG VACCINE (JAPAN)
FOR INTRADERMAL USE

BK

PANTONE 301 CVC

PANTONE Warm Gray 5 CVC