

Duosol with 2 mmol/l Potassium Solution for haemofiltration

2 x 5000 ml (two-chamber bags, 4445 ml and 555 ml) | Art. 4509 | Intravenous use 07 4509/0121

Warning: Use only after opening the peel seam and mixing of the two solutions.

Before mixing, the bicarbonate solution (large chamber) contains (4445 ml):	Do not store above 25 °C. Do not store below 4 °C.	1000 ml ready-to-use solution for haemofiltration contain [mmol/l]:
Sodium chloride 27.47 g	Free from bacterial endotoxins.	Na ⁺ 140
Sodium hydrogen carbonate 15.96 g	Use only if the container, peel seam and connectors are undamaged and intact and if the solution is clear and colourless and free from visible particles. Remove outer wrap only immediately before use.	K ⁺ 2.0
Carbon dioxide, water for injections	Read the leaflet for the shelf life after mixing.	Ca ⁺⁺ 1.5
	Discard any unused contents. Read the package leaflet before use.	Mg ⁺⁺ 0.5
	Keep out of the sight and reach of children.	Cl ⁻ 111
Before mixing, the electrolyte solution (small chamber) contains (555 ml):	Jauhkan dari kanak-kanak.	HCO ₃ ⁻ 35.0
Sodium chloride 2.34 g	Controlled Medicine.	Glucose anhydrous 5.6 (equiv. to 1.0 g)
Potassium chloride 0.74 g	Marketing Authorisation number: MAL09011488AZ	Theoret. osmolality [mOsm/l] 296
Calcium chloride dihydrate 1.10 g		pH 7.0-8.0
Magnesium chloride hexahydrate 0.51 g		
Glucose monohydrate 5.49 g		
equiv. 5.0 g glucose anhydrous		
Hydrochloric acid 25% (for pH adjustment), water for injections		

POM

Manuf. date:

Batch no.:

Expiry date:

GTIN: 04046963670283

Approval for Printing

B | BRAUN Avitum AG

Approved for Printing



Date Signature

Name in capital letters

Manuf. date: DD.MM.YYYY

Batch no.: xxxxxxxx

Expiry date: DD.MM.YYYY

The Manufacturing Date, the Batch Number, the Expiry date and a 2D Datamatrix-Code (carrying the Batch Number, Expiry Date and the GTIN) will be printed in the green marked blank field during production.

B. Braun Avitum AG
34212 Melsungen
Germany

Manufacturing site:
B. Braun Avitum AG
Kattenvenner Str. 32
49219 Glandorf, Germany

Product registration holder:
B. Braun Medical Industries Sdn. Bhd.
11900 Bayan Lepas, Penang
Malaysia

B | BRAUN