



Nutriflex[®] plus


Controlled Medicine

Intravenous Infusion Solution for Parenteral Nutrition

The mixed solution contains:

Total Amino Acids 48.1 g/l
Anhydrous Glucose 150 g/l

Nitrogen 6.8 g/l
Osmolarity 1400 mOsm/l
kJ (kcal)/l 3310 (790)



Sterile and apyrogenic.
Do not use unless solution is clear and container is intact. Contents of the two chambers must be mixed just before use. Discard any unused solution. Read instruction leaflet before use. Store below 30 °C, protected from light and out of the reach of children.

1. Unfold the bag and place it on a sturdy surface.
2. Open the peel seam by pressing with both hands on one chamber of the bag.
3. The homogeneous mixture is ready for use.

1500

1250

Manuf. Date

Batch No

Exp. Date

1000

750

500

1000 ml of the ready-to-use infusion solution contains after mixture of both chambers:

Upper chamber contains (400 ml):

Isoleucine 2.82 g
Leucine 3.76 g
Lysine Hydrochloride 3.41 g
△ Lysine 2.73 g
Methionine 2.35 g
Phenylalanine 4.21 g
Threonine 2.18 g
Tryptophan 0.68 g
Valine 3.12 g
Arginine Monoglutamate 5.98 g
△ Arginine 3.24 g and Glutamic Acid 2.74 g
Histidine Hydrochloride Monohydrate 2.03 g
△ Histidine 1.50 g
Alanine 5.82 g
Aspartic Acid 1.80 g
Glutamic Acid 1.47 g
Glycine 1.98 g
Proline 4.08 g
Serine 3.60 g
Magnesium Acetate Tetrahydrate 1.23 g
Sodium Acetate Trihydrate 1.56 g
Sodium Dihydrogen Phosphate Dihydrate 3.12 g
Potassium Hydroxide 1.40 g
Sodium Hydroxide 0.23 g

Water for Injections 360.2 g

Excipient: Citric Acid

Lower chamber contains (600 ml):

Glucose Monohydrate 165 g
△ Glucose 150 g
Calcium Chloride Dihydrate 0.53 g
Water for Injections 489.4 g
Excipient: Citric Acid

1000 ml contains:

Electrolytes:

Na ⁺	37.2 mmol/l	Cl ⁻	35.5 mmol/l
K ⁺	25.0 mmol/l	H ₂ PO ₄ ⁻	20.0 mmol/l
Ca ⁺⁺	3.6 mmol/l	Acetate ⁻	22.9 mmol/l
Mg ⁺⁺	5.7 mmol/l		

Jauhkan dari kanak-kanak

Content: **1000 ml**

US Patent Number: 4.917.925

Reg. No.: MAL19962093AZ

Product Registration

Holder:

B. Braun Medical
Industries Sdn. Bhd.
Bayan Lepas Free
Industrial Zone
11900 Bayan Lepas
Penang, Malaysia

Manufactured by:

B. Braun Medical SA
1023 Crissier
Switzerland



Infusion port



Additive port



715/CH15CR1378/0416

B | BRAUN

HKS 52 Schwarz
148x280 mm
715/CH15CR1378/0416

Malaysia
Font size 8

