

Nutriflex® Lipid Peri Novo

0846/12631303/1124

Emulsion for Infusion
1250 ml

MFG date
Lot:
EXP

955 kcal	Lipid 50 g
Amino acid 40 g	Carbohydrate 80 g



Excipients
Citric acid monohydrate (for pH adjustment)
Glycerol
Egg phospholipids for injection
Sodium oleate
all-rac-alpha-Tocopherol
Water for injections

Electrolytes:

Sodium	50 mmol	Energy in the form of lipids	1 990 kJ (475 kcal)
Potassium	30 mmol	Energy in the form of carbohydrates	1 340 kJ (320 kcal)
Magnesium	3.0 mmol	Energy in the form of amino acids	670 kJ (160 kcal)
Calcium	3.0 mmol	Non-protein energy	3 330 kJ (795 kcal)
Zinc	0.03 mmol	Total energy	4 000 kJ (955 kcal)
Chloride	48 mmol		
Acetate	40 mmol		
Phosphate	7.5 mmol		

Amino acid content:	40 g	Osmolality	950 mOsm/kg
Nitrogen content:	5.7 g	Theoretical osmolality	840 mOsm/l
Carbohydrate content:	80 g	pH	5.0 – 6.0
Lipid content:	50 g		

0/12613265/0417

1250 ml of the ready-for-use emulsion contain after mixing of the three chamber contents from this chamber (250 ml):

Soya-bean oil, refined	25.00 g
Medium-chain triglycerides	25.00 g

Emulsion for infusion
Intravenous use
Read the package leaflet before use. Keep out of the sight and reach of children. Jauhi ubat dari kanak-kanak.
High in sodium – see leaflet for further details. Contain egg phospholipid from hens. After removing the protective overwrap and after mixing of contents of the bag (Before spiking of the infusion port). Chemical and physicochemical in-use stability of the mixture of amino acids, glucose and fat was demonstrated for 7 days at 2-8 °C and additional 2 days at 25 °C. After admixture of compatible additives. From a microbiological point of view, the product should be used immediately after admixture of additives. If not used immediately after admixture of additives, in-use storage times and conditions prior to use are the responsibility of the user. After first

peripheral
or central
venous


MCT/LCT

with
Electro-
lytes

1000 mL



opening (spiking of the infusion port). The emulsion is to be used immediately after opening of the container.
Do not store above 25 °C. Do not freeze. If accidentally frozen, discard the bag. Keep the bag in the protective overwrap in order to protect from light.
Reg Number: MAL24116026AZ
Controlled Medicine
Ubat Terkawal
Prescription Only Medicine



1250 ml of the ready-for-use emulsion contain after mixing of the three chamber contents from this chamber (500 ml):

Glucose monohydrate equivalent to glucose	88.00 g
Sodium dihydrogen phosphate dihydrate	1.170 g
Zinc acetate dihydrate	6.600 mg



250 mL


1250 ml of the ready-for-use emulsion contain after mixing of the three chamber contents from this chamber (500 ml):

Isoleucine	2.340 g
Leucine	3.130 g
Lysine hydrochloride equivalent to lysine	2.840 g
Methionine	2.273 g
1.960 g	


Phenylalanine	3.510 g	Glutamic acid	3.500 g
Threonine	1.820 g	Glycine	1.650 g
Tryptophan	0.570 g	Proline	3.400 g
Valine	2.600 g	Serine	3.000 g
Arginine	2.700 g	Sodium hydroxide	0.800 g
Histidine hydrochloride monohydrate	1.690 g	Sodium chloride	1.081 g
equivalent to histidine	1.251 g	Sodium acetate trihydrate	0.544 g
Alanine	4.850 g	Potassium acetate	2.943 g
Aspartic acid	1.500 g	Magnesium acetate tetrahydrate	0.644 g
		Calcium chloride dihydrate	0.441 g

B|BRAUN

Product registration holder:
B. Braun Medical Industries Sdn. Bhd.,
Bayan Lepas Free Industrial Zone,
11900 Bayan Lepas, Penang, Malaysia.



Infusion port



Additive port

Product Owner:
B. Braun Melsungen AG,
Carl-Braun-Straße 1,
34212 Melsungen, Germany.

Manufactured by:
B. Braun Melsungen AG,
Am Schwerzelshof 1,
34212 Melsungen, Germany.

Approval for Printing

B|BRAUN Melsungen AG

Approved for Printing

Approved for Printing when corrected

New draft required

LRA:

Date: _____ Signature: _____

Name in capital letters _____

LMS:

Date: _____ Signature: _____

Name in capital letters _____

Black
HKS 46 (colour cliche)

Bag outline (no printing color! For information only!)

Dimension (printable area): 250 x 265 mm
LLD-Spec.: M53
Serialized: _____

LIFE_N__3CB__MY__0846
Construct: 0846/12631304/1124
BLACK (cliche): 0846/12631303/1124
Colour (cliche): 0/12613265/0417

Product: LIFE N 3-Chamber bag 1250 ml
Production site: Melsungen (LIFE N)

Font size: 8,0 pt.

V-0013

File name:
12631304_NutriflexLipidPeriNovo_LIFEN--3CB-1250ml__MY.indd

Date of last saving:
13.11.2024 - 11:45:39

Example for variable online print with dummy content:

? n.k.