

GTIN: XXXXXXXXXXXXX
S.N.: XXXXXXXXXXXXX
LOT: XXXXXXXXX
MFD: MM.YYYY
EXP: MM.YYYY
Batch details along with variable data shall be over coded with 2D during online packing. Whereas 2D code contains GTIN, S.N, Lot Number & Expiry date



MSN LABORATORIES PRIVATE LIMITED		PACKAGING DEVELOPMENT	
Artwork Information		Specification for Printed Carton	
		Test	Specification
Brand Name	LURATA 80	Substrate	Folding Box Board
Generic Name	Lurasidone Hydrochloride Tablets 80 mg	GSM	270 ± 5% GSM
Pack Style / Counts	1 x 10's	Varnish/Lamination	UV varnish except over printing zone
Auto Cartonator / Manual	Auto Cartonator	Mode of supply	Bundles 25's, 50's & 100's
Reverse Tuck / Lock Bottom	Reverse Tuck	Pharma Code	2801
Dimension	(L)122 x (W)62 x (H)25 ± 1 mm	Special Features	NA
Item Code	D022777-00	Font Type	Arial Narrow, OCR-B 10 BT
Supersede Code	NA	Font Size	6 Pt
Version	01	Developed by	G. Mahesh
Date & Time	08.01.2024 & 15:30	Reviewed by	B. Veeranna
Country	Malaysia	Colours	P 300 C & 30% of P 300 C P 485 C P 273 C Black Outer lines not to be printed (Unvarnished area) Outer lines not to be printed
Customer	Mansa Healthcare		
Instruction: Carton should be supplied suitable to run on auto cartonator as per the auto cartonator drawing attached Overlapping of Closing/Collar flaps of carton to be avoided.			

during online packing,
Whereas 2D code contains GTIN,
S.N, Lot Number & Expiry date

MSN

LURATA 80
Lurasidone Hydrochloride 80 mg
Film-coated Tablets

FOR ORAL USE
30 tablets

Ubat Terkawal/Controlled Medicine

MSN

LURATA 80

0022778-00

Unvarnished Area
30 x 62 mm

LURATA 80

Each film coated tablet contains
Lurasidone hydrochloride 80 mg
Product contains Lactose

DOSAGE: As directed by the Physician

Distribution Category: POM

Do not store above 30°C.
Protect from light and moisture.

Jauhi dari kanak-kanak
Keep out of the sight and reach of children

Please Refer Packaging Insert
for detailed information.

Product Registration Holder
Mansa Healthcare Sdn Bhd
B-05-3A, 3 Two Square, No 2, Jalan 19/1, 46300
Petaling Jaya, Selangor

Manufactured by:
MSN Laboratories Private Limited,
Formulation Division, Unit-II,
Sy.No. 1277, 1319 to 1324,
Nandigama (Village & Mandal),
Rangareddy (District),
Telangana - 509 228, India.

M.L.No.: 5/MN/TS/2014/F/G
MAL No.: MAL23106004AZ

8 904159 684000

MSN


LURATA 80

30 mm 122 mm 30 mm

62 mm 30 mm 62 mm 30 mm

MSN LABORATORIES PRIVATE LIMITED		PACKAGING DEVELOPMENT	
Artwork Information		Specification for Printed Carton	
		Test	Specification
Brand Name	LURATA 80	Substrate	Folding Box Board
Generic Name	Lurasidone Hydrochloride Tablets 80 mg	GSM	270 ± 5% GSM
Pack Style / Counts	3 x 10's	Varnish/Lamination	UV varnish except over printing zone
Auto Cartonator / Manual	Auto Cartonator	Mode of supply	Bundles 25's, 50's & 100's
Reverse Tuck / Lock Bottom	Reverse Tuck	Pharma Code	2802
Dimension	(L)122 x (W)62 x (H)30 ± 1 mm	Special Features	NA
Item Code	D022778-00	Font Type	Arial Narrow, OCR-B 10 BT
Supersede Code	NA	Font Size	6 Pt
Version	01	Developed by	G. Mahesh
Date & Time	08.01.2024 & 15:20	Reviewed by	B. Veeranna
Country	Malaysia	Colours	
Customer	Mansa Healthcare		

Instruction: Carton should be supplied suitable to run on auto cartonator as per the auto cartonator drawing attached
 Overlapping of Closing/Collar flaps of carton to be avoided.



MSN LABORATORIES PRIVATE LIMITED		PACKAGING DEVELOPMENT	
Artwork Information		Specification for Printed Carton	
Brand Name	LURATA 80	Test	Specification
Generic Name	Lurasidone Hydrochloride Tablets 80 mg	Substrate	Folding Box Board
Pack Style / Counts	10 x 10's	GSM	300 ± 5% GSM
Auto Cartinator / Manual	Auto Cartinator	Varnish/Lamination	UV varnish except over printing zone
Reverse Tuck / Lock Bottom	Reverse Tuck	Mode of supply	Bundles 25's, 50's & 100's
Dimension	(L)122 x (W)62 x (H)65 ± 1 mm	Pharma Code	2803
Item Code	D022779-00	Special Features	NA
Supersede Code	NA	Font Type	Arial Narrow, OCR-B 10 BT
Version	01	Font Size	6 Pt
Date & Time	08.01.2024 & 15:10	Developed by	G. Mahesh
Country	Malaysia	Reviewed by	B. Veeranna
Customer	Mansa Healthcare	Colours	P 300 C & 30% of P 300 C P 485 C P 271 C Black Order lines not to be printed
Instruction: Carton should be supplied suitable to run on auto cartinator as per the auto cartinator drawing attached. Overlapping of Closing/Collar flaps of carton to be avoided.			